



Features

1. Luminance: ≥ 140000 Lux
2. Color temperature range: 4500 ± 500 K
3. Spot diameter: 150 ~ 300 mm
4. Illumination depth: ≥ 900 mm
5. Total irradiance: ≤ 1000 W/m
6. Temperature at surgeon's head: $\leq 1^\circ$
7. Power supply: ~ 220 V 50 Hz
8. Input power: 200 VA
9. Color indication index: $\geq 97\%$
10. Bulb life: ≥ 50000 h
11. Six-dimensional joint: 360°
12. Lowest height of installation: 2600 mm