

ICU Infant Care

TRBI 100[®] Infant Incubator

Environmental temperature

Operating range: 20 °~30°

(Set temperature should be higher 3° than the environment temperature)

Ambient air movement rate: <0.3m/s

Specification

Power supply: AC220V-230V/50Hz

Max Power Output: 850VA

Mode of control: Air mode control and skin mode control by micro-computer

Air temperature control range: 25°-37°

Skin temperature control range 34° - 37°

Temperature variability: ≤0.5°

Temperature uniformity: ≤0.8°

Noise inside hood: ≤60dB(A) (Environment noise≤45dB(A))

Failure alarm: over temperature alarm, deviation alarm, sensor failure alarm, fan motor alarm, power failure alarm, system alarm and so on.

