

ANY CARE



About DS maref

DS MAREF is the leading company of various Pneumatic medical systems and exporting our products to 73 countries since 1986.

Through our main marketing fields such as Medical, Esthetic, Obesity, Sports/Rehabilitation and Home care, we have supplying various professional equipments to each field. Besides, we are participating in international medical, beauty and health care exhibitions by acquiring the certificates of FDA, CE, ISO9001, KGMP, KFDA, SFDA and etc. Even in competitiveness of the products, we have securing unequalled competitiveness rather than others in quality, design and price.

We will listen to more customers' ideas, and finally, do the best to be a top company of No.1 Pneumatic Compression Therapy Systems in the world through constant developing new products.

ISO 13485

CE

FDA

SFDA

KFDA

KGMP

IPC Therapy system

ANVISA

TGA

INMETRO

Certificates & Awards

2011 President Award
2010 Anvisa in Brazil
2007 Authentication of Merit certificate
2006 Innovation management awards
2006 KOTRA B2B e-Trade awards
2006 KFDA awards by Prime minister
2006 KGMP registration
2004 The 34 precious Technique bronze awards
2004 Korean World class products award
2004 Vice president award
2004 SFDA in China
2003 Success Design product
2003 FDA in USA
2002 ISO9001/EN13485
2002 CE marks (DL / MK series, DL1200, DVT-2600)
1986 DS MAREF is established

Our activities

Overseas activities

- Above 15 International Exhibitions
- Exporting to 73 Countries (2011)
- Sales increases rate, 30~40%
- FDA, CE, ISO13485, SFDA, TGA, KGMP



2011 독일 메디카전시

Domestic activities



- Market share No.1 in medical, in esthetic, home care fields
- Advertisement in KIMES show, Online shopping malls, TV, Radio, Commercial Catalogue



2012 키메스 의료전시

DoctorLife air compression system

- All human bodies are managed by smooth circulation of blood and lymph.
- DOCTOR LIFE drainage system system helps that various patients are recovered from their disease by promoting blood and lymph circulation and metabolism.



DOCTOR LIFE *ANYCARE*



ANYCARE offers easy operation and portable type with more polished design as a home care device. This design came from the concept of women's Luxury bags.

V6_pink & white

ANYCARE

DOCTOR LIFE®

TOTAL SOLUTION BRAND SINCE 1986

TY




ssure
4

controller 2

controller 1

ANYCARE



ANCYARE premium features

ANYCARE premium Specifications	
Voltage	AC220V
Hz	50/60Hz
Power consumption	32W(40VA)+10%
Pressure range	Max. 240mmHg±15%
Time range	0~30분
Pressure control	Individual pressure setting




ANYCARE Specifications	
Voltage	AC220V
Hz	50/60Hz
Power consumption	32W(40VA)+10%
Pressure range	Max. 240mmHg±15%
Time range	0~30 minutes

Cuff features

ANYCARE basic set includes 1 machine, 2 leg cuffs, 1 double hose, 1 single hose, 2 foot pad

Options: Arm cuff, Waist cuff, Center body cuff, Pelvis cuff, Leg extension zipper, Arm extension zipper



Thank you



DS MAREF
주식회사 대성마리프

(주)대성마리프 | 경기도 군포시 공단로 298-24
Tel : (031) 459-7211 Fax : (031) 459-7215
E-mail : ds@dsmaref.com www.dsmaref.com