

Beurer BC 28 wrist blood pressure monitor

Straightforward and easy to use



Simple and easy to use – that is what makes our wrist blood pressure monitors stand out. These handy devices prove themselves to be particularly reliable companions on the move, meaning you can monitor your blood pressure regularly.

- Fully automatic blood pressure and pulse measurement on the wrist
- Average of all saved measured values
- Average of morning and evening blood pressure for the last 7 days
- Clear display
- 2 x 60 memory spaces
- Risk indicator
- Arrhythmia detection
- Medical device
- For wrist circumferences from 14.0 to 19.5 cm
- Date and time/automatic switch-off
- Incorrect usage message
- Low battery indicator
- Storage box

Product designation	Wrist blood pressure monitor
Automatic switch-off	yes
Number of memory spaces	2 x 60
Arrhythmia detection	yes
Storage bag/box included in delivery	yes
Batteries	yes

Low battery indicator	yes
Date and time	yes
Calculating the average value	yes, 7 days
Large, easy-to-read display	yes
Cuff size in cm	14.0 – 19.5
Medical device	yes
Incorrect usage message	yes
Risk indicator	yes
Product dimensions (L x W x H)	84 x 60 x 29 mm