

Beurer BF 105 Body Complete diagnostic bathroom scale

Detailed full-body analysis



The diagnostic bathroom scale is a true all-rounder. Use the scale for a complete analysis of your body values. The graphical evaluation and trend display of your data helps you to achieve your personal goals.

- Detailed full-body analysis
- Precise upper/lower body measurement with 8 electrodes
- Removable remote handheld unit with blue illumination
- Weight, body fat, body water, muscle percentage, bone mass, AMR/BMR calorie display
- Precise upper/lower body measurement
- Removable XXL display: 120 x 46 mm
- State-of-the-art dot matrix display
- Intuitive user menu in: D, GB, F, ES, I
- Large-surface electrodes made of brushed stainless steel
- USB and *Bluetooth*® transfer to PC
- User memory spaces: 10 + guest mode
- Measurement memory spaces: 30
- Switch-on technology: Button-on
- Automatic switch-off, overload indicator
- Weight capacity: 180 kg
- Digit size: 22 mm

- "beurer HealthManager" software and app available for free download



User memory spaces

10 user memory spaces

Weight capacity

Up to 180 kg



Body analysis

Weight, body fat, body water, muscle percentage, bone mass; calorie display

Full-body analysis

Precise upper and lower body measurement

AMR/BMR

Calculation of BMR (Basic Metabolic Rate), which denotes the basic level of metabolism, and AMR (Active Metabolic Rate), which denotes the total metabolic rate, including exercise and physical activity



HealthManager

Compatible with the "beurer Health Manager" app

Modern, networked health management! Beurer's specially developed 'beurer Health Manager' app enables you to keep an eye on your data at all times – from blood pressure and blood glucose to activity, heart rate, oxygen saturation and weight.

Trend display

Graphical evaluation and trend display

XXL format

All measurements: Weight, body fat, body water, muscle mass, bone mass is shown simultaneously in the unique display

Product designation	Diagnostic bathroom scale
Connect	yes
5 activity levels	yes
Batteries	6 x 1.5 V AA batteries
Illuminated display	yes
User memory spaces	10
Beurer Connect	yes
Electrodes	in brushed stainless steel
Standing surface size in cm	33 x 30
AMR calories	yes
BMR calories	yes
Bone mass	yes
Compatible from	iOs 10.0 and Android 5.0 by Bluetooth® 4.0
Body fat display	yes
Body water display	yes
Measurement memory spaces	30
Muscle percentage display	yes
Product dimensions in cm	33 x 33 x 4
Weight capacity/graduation	180 kg / 100 g
Trend display	in app
Transfer with	USB and Bluetooth®
Conversion	kg/lb/st
Digit size in mm	22