

# Beurer BG 13 diagnostic bathroom scale



Weighing made easy! The scale with large LCD display shows you your weight, body fat, body water and muscle percentage and even calculates your BMI!

- With large LCD display
- Weight, body fat, body water, muscle percentage
- With BMI calculation
- Subtle striped design
- XL display
- Electrodes made of brushed stainless steel
- 5 activity levels
- User memory spaces: 10
- Switch-on technology: Quick Start
- Automatic switch-off, overload indicator
- Digit size: 38 mm



## **User memory spaces**

10 user memory spaces

## **Weight capacity**

Up to 150 kg



## **Body analysis**

Weight, body fat, body water, muscle percentage, bone mass; calorie display

## **BMI**

The scale calculates your Body Mass Index and shows it on the display

## **Slim design**

Slim design only 18 mm

## **XL display**

The product features a convenient XL display

Product designation	Diagnostic bathroom scale
5 activity levels	yes
Batteries	1 x 3 V CR2032 battery
User memory spaces	10
BMI calculation	yes
Switch-on technology	Quick Start
Electrodes	in brushed stainless steel
Standing surface size in cm	30 x 30
Body fat display	yes
Body water display	yes
Muscle percentage display	yes
Product dimensions in cm	30 x 30 x 1.8
Weight capacity/graduation	150 kg / 100 g
Standing surface made of safety glass	yes
Conversion	kg/lb/st
Digit size in mm	38