

Beurer FT 65 multi functional thermometer

Suitable for the whole family



A real all-rounder! In addition to ear and forehead measurement, the thermometer can also determine surface temperatures. The temperature alarm and the date and time display are perfect additions for simple and convenient measurement.

- 6-in-1 function: ear, forehead and surface temperature, temperature alarm, date and time, 10 memory spaces
- Up to 37.9°C = green, from 38.0°C = red
- Measurement in seconds
- Infrared measurement
- LED temperature alarm from 38°C
- Displays measurements in °C and °F
- No mercury, without glass
- Can be disinfected
- 10 memory spaces, with date and time
- Automatic switch-off
- Medical device
- Large display

Automatic switch-off	yes
Low battery indicator	yes
Date and time	yes

Medical device	yes
Measurement of the surface/object/liquid temperature	yes
Measurement of ear temperature	yes
Measurement of forehead temperature	yes
Visual temperature alarm (coloured LED)	yes
Acoustic signal	yes
Memory spaces	10
Can be set to °C/°F	yes
Product weight	approx. 90 g
Product dimensions	38.2 x 138 x 46.5 mm