

# Beurer MG 55 tapping massager

## Penetrating and relaxing



With the penetrating and relaxing tapping massage you can relax and loosen your muscles.

- Penetrating and relaxing tapping massage
- For relaxing and loosening muscles
- Continuously adjustable massage intensity
- Heat function can be selected separately
- 3 interchangeable massage attachments (nub attachment, round attachment, cylinder-shaped attachment) with threads for secure hold and easy use
- Non-slip handle
- Ergonomically shaped

Product designation	Tapping massager
Function levels	1
Output in watts	20
Type of massage	Tapping massage
Product dimensions	9,5 x 43,5 x 9,5 cm