

# Beurer PO 30 pulse oximeter

**Small and handy**



You can determine the oxygen saturation in your blood with the pulse oximeter. This may be necessary during sport at high altitudes or if you have health problems. Your pulse is also shown graphically on the colour display.

- Measurement of arterial oxygen saturation (SpO<sub>2</sub>) and heart rate (pulse)
- Particularly suitable for persons with: Heart failure, Chronic obstructive pulmonary diseases, Bronchial asthma
- Also suitable for sports at high altitudes
- Colour display with 4 available views
- Adjustable display brightness
- Graphic pulse display
- Small and light – for use at home and on the move
- Automatic switch-off
- Low battery indicator
- Incl. strap and belt pouch

Product designation	Pulse oximeter
Available views	4
Battery-operated	yes
Batteries	2 x 1.5 V AAA batteries
Colour display	yes
Graphic pulse display	yes

Medical device	yes
Product weight	approx. 57 g
Product dimensions	61 x 36 x 32 mm