

Five Function Electric Bed KS-858

Kcare



This electric five Function hospital Bed is carefully designed for patients during their treatment. The head and foot board can be removed for the convenience of patients' transfer. The side rails are made of ABS material,

can be easily rise and lower. With 4 pcs motors under the bed. With 4 pcs castors with brakes, this is convenient for transportation.

Product Description

Product Name	Five function electric bed
Model	GS-858
Dimension:	2260*1060*(430~720) mm<L*W*H>
Motor	Linak motor
Bed board	1.Stainless Steel material, 2. Overload capacity: 175KG.Easy to clean.
Guardrail	1.PP material, 2. Interspace figure: a. Between the Guardrails: less than 3cm, b. Between the Guardrail and Head/Foot board: less than 6cm.
Head& Foot board	1.ABS material, 2.7KG (15lb), 3. Smooth and easy to clean.
	10 segments sleeping board,5" central lock castor, back CPR manual, IV pole
	ABS side rail with function button. LINAK motor with battery
Castor	5 inches castor, Secured Brand.
Functions	1.Back section up& down, 2. Knee section up& down, 3. The whole bed up& down, 4. Trendelenburg and reverse Trendelenburg, 5. Back section decompression function.

Standard Accessory

Item Name	Quantity
Head and Foot Board	1 pair.
Guardrail	4 pieces.
Motor	4 sets.

Remote	1 set.
Angle Indicator	4 sets.
Castor	4 pieces.
Side brake pedal	2 pieces.
I.V drip hole	4 sets.
Drainage bag hook	4 pieces.
CPR	1
UPS (back-up battery)	1