



TRCB 3004 Multi-purpose obstetric table

TRCB 3004 multi-purpose obstetric table is designed for delivery, abortion and gynecological examination. It is characterized by easy operation, safe and reliable, economical use. It is composed of bed surface, bedstead and bed

base. The surface includes back section, seat section and leg section. The up-down movement of back section and forward and back tilt if bed surface are controlled by the operating hand wheel, which makes doctors and patients feel comfortable.



Remark

- 1.leg table with drawer
- 2.back section, forward and backward are adjusted by the machine
- 3.up and down of the table movement by hydraulic
- 4.the cover of oil pump and bed base are made of stainless steel

Main specification parameters:

Length and width of table	1800mm*600mm
Height adjustment range	Min 650mm* Max 895mm
Size of sections
Back section	600mm*800mm
Seat section	600m*400mm
Leg section	540mm*500mm
Max. degree of back rest adjustment	up>75°down>6°
Max. degree of forward and backward tilting of tabletop	forward >10° backward >25°