

## ZT-200 Weight & height measuring scale



### Specification

Floor mounted for adult

Weight scale must be placed in front of person at standard height

Mechanical type

The Weighing scale is equipped with pulley which can make the movement in short distance more convenient

Capacity: 200 kg

Scale deviation : 500gm

Powder Coated Painting

Unit: kg & lb and meter & feet for weight & height respectively

Height measurement range : 70-190cm height.

Height per division : Max. 0.5 cm

Packing : Export Standard Packing

Net weight 11 kg (approx.)

Gross weight 12.5 kg (approx.)