



ITEM	SPECIFICATION
WEIGHT CAPACITY	350lbs (160kgs)
DRIVE WHEEL	13" x 4" (330mm x 100mm)
FRONT CASTER(WHEEL)	8" x 2" (200mm x 50mm)
REAR CASTER(ANTI-TIPPER)	6" x 1.375" (155mm x 35mm)
MAX SPEED	5MPH (8KPH)
BATTERY SPECIFICATIONS	12V 50Ah x 2pcs
BATTERY RANGE	25miles (40km)
CHARGER TYPE	5 Amp Off Board
CONTROLLER TYPE	VSI 50Amp
MOTOR TYPE	450W x 2pcs
LENGTH	44" (1110mm)
WIDE	24.5" (620mm)
WEIGHT:W/BATTERY	213lbs (97kgs)
SUSPENSION	Rear Caster
TURNING RADIUS	26" (60mm)
HEIGHT	45" (1140mm)
SEAT WIDTH	20" (508mm)
SEAT HEIGHT	21" (535mm)
SEAT DEPTH	18" (460mm)
BACK HEIGHT	27" (685mm)
WHEEL BASE	12.5" (320mm)
LEG ROOM	16.5" (420mm)