



MIKO is a piston compressor nebulizer for aerosoltherapy. Designed for home-care with compact size and carrying bag for easy transport of unit, accessories and medications.

#### STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
- Adult Masks with elastic strings
- Child Masks with elastic strings
- Mouthpiece
- Nosepiece (non invasive)
- Air Tube
- Air Filter
- Carrying Bag (MIKO model only)

#### TECHNICAL FEATURES

Compressor	Oilless and maintenance-free piston pump		
Power Feeding	230V-50Hz - other configurations on request		
Power consumption	170 VA		
Max Pressure	2.2 bar	32 psi	220 kPa
Operating Pressure	0,95 bar	13.8 psi	95 kPa
Max Air Flow	14 l/min		
Operating Air Flow	4.8 l/min		
Neb-rate with 4 ml solution	0,35 ml/min		
Mass Median Aerodynamic Diameter	2,44 $\mu$ m (MMAD) according to EN-13544-1		
Respirable Fraction (< 5 $\mu$ m)	> 80%		
Duty Cycle	Non-stop operation		
Noise Level (EN 13544-1 Standard)	Approx. 56 dBA		
Weight	1,45 kg		
Size (Main Unit)	cm 14 x 17 x 11		
Years of Warranty	3		
Shipping carton	6		
Place of Manufacturing	Italy		

3 YEARS  
WARRANTY

Italian  
Medical  
Touch  
MADE IN ITALY

IP21

#### CODE

RE 300600/03	MIKO- with carrying bag
RE 300600/12	MIKO BASIC - without carrying bag

CE 0123