

Two crank Manual Bed ks-332

Kcare



This two Function Manual Bed is carefully designed for patients during their treatment. The head and foot board can be removed for the convenience of patients' transfer. The side rails are made of aluminum alloy. And the gap between two columns is according to the international standard. All the position is achieved by the foldable and concealed manual cranks. This will

provide long use life of cranks. With 4 pcs castors with brakes, this is convenient for transportation.

Product description:

Description of product:

1. Give the patient comfortable, warm and home felt by using the ABS head and foot board.
2. Button type Aluminum guard rail with 6pieces of stainless-steel bracket. Built-in control button, easy to control.
3. Deluxe 5 inches double face casters

Product parameters

Name	manual two crank care bed
Item NO.	ks-332
size	L2120*W970*H485mm
standard accessories	impendent brake castors
	I.V. stand*1; I.V. stand hole×4 drainage hook
products functions	backrest: $75^{\circ}\pm 5^{\circ}$; footrest: $45^{\circ}\pm 5^{\circ}$
packing size	2040*940*230
Bed panel	10 segments sleeping board Stainless steel
Raw material	steel
Max load	250kgs