

## **One Crank Manual Bed KS-S202yh**

### **Kcare**

This one Function Manual Bed is carefully designed for patients during their treatment. The head and foot board can be removed for the convenience of patients' transfer. The side rails are made of aluminum alloy. And the gap between two columns is according to the international standard. All the position is achieved by the foldable and concealed manual cranks. This will provide long use life of cranks. With 4 pcs castors with brakes, this is convenient for transportation.





## **Product description:**

Description of product:

1. Give the patient comfortable, warm and home felt by using the ABS head and foot board.
2. Button type Aluminum guard rail with 6pieces of stainless-steel bracket. Built-in control button, easy to control.

### 3. Deluxe 5 inches double face casters

#### Product parameters

Brand Name :	K-care
Product Name :	one crank manual Hospital Bed
Place of Origin :	Guangdong, China
Model Number:	KS-S202yh
Size:	L2150*W980*H500mm
Function:	Backrest 0~75°±5°
Head and foot board :	1 pair, ABS material, with patient name card
Color :	Blue/brown
Wheels:	5", 4pcs, with individual brake
Guardrail	1pair, aluminum
I.V. Stand	1pc, S.S, telescopic, with 2 hooks
I.V. Stand hole	6pcs
Drainage hook	4pcs
Bed panel	10 segments sleeping board Stainless steel
Packing sizes:	2050*960*230mm
Assemble :	Frame has assembled, Side rails/ head and
Package:	Carefully selected high-standard packaging one carton