

Three Function Electric Bed KS-888

Kcare





Powered by DIYTrade.com

This Three Function Electric Bed is carefully designed for patients during

their treatment. The head and foot board can be removed for the convenience of patients' transfer. The side rails are made by Aluminum side rail. And the gap between two side columns is according to the international standard. All the position is controlled by electric, use a remote to control all the function. This will provide long use life of cranks.

With 4 pcs castors with individual brakes, this is convenient for transportation.

Product parameters

Model	KS-888
Side	2120x970x210--600mm
Bed panel	Stainless steel
Head and foot board	ABS
Side rail	Aluminum alloy
Wheel	Central Lock Brake
Motor	Linak with battery
	head side handling bar
Bed panel	10 segments sleeping board Stainless steel
Backrest function	75°±5°
Footrest function	45°±5°
The whole bed up & down	430/720mm
Maximum load	250kgs no movement – 175kgs in movement

Standard device

ABS head and foot board	1 pair
Side rail	1 pair
Motor	3 pcs
5-inch wheel	4 pcs
pedal	1 pc
IV Stand	1 pc

IV drip holes	4 pcs
Drainage bag hooks	4 pcs