

# NIDEK *Medical*

**Proudly Presents  
The Most Powerful  
Member of our  
Nuvo Family of  
Oxygen Concentrators**



.....Mark 5  
**Nuvo 8**

- Liter Flow -- 2 to 8 LPM
- Ideal for High Flow Patients
- Low Maintenance Costs
- Made in the USA

**CE**  
0413





# NIDEK

*Medical*

## Additional Features:

- Built on the Reliable Nuvo Platform
- Easy Access to filter, battery, and hourmeter
- Superior Grade of Molecular Sieve
- Spring Mounted High Efficiency Thomas Nexus Compressor

## Specifications\*

**Nuvo 8 Oxygen Concentrator**

**Nuvo 8 Oxygen Concentrator  
with Oxygen Concentration Status Indicator (OCSI)**

**115 Volt - 60 Hz**

**580**

**585**

**230 Volt - 50/60 Hz**

**980**

**985**

Average Power Consumption  
Oxygen Concentration (at 8 LPM)  
Liter Flow  
Outlet Pressure  
Weight  
Dimensions

500 Watts  
90% (+5% / -3%)  
2 to 8 LPM  
17 PSIG (1.2 BAR)  
53.5 lbs / 24.3 kg  
15.5"W x 15.6"D x 27.8"H

490/625 Watts  
90% (+5% / -3%)  
2 to 8 LPM  
17 PSIG (1.2 BAR)  
55.5 lbs / 25.2 kg  
39.4 cm x 39.6 cm x 70.6 cm

## Accessories and Re-order Items

Compressor Intake Filter  
External Cabinet Filter  
Humidifier  
Cannula with 7' tubing  
Nuvo 8 User's Guide

9250-1180  
9250-1025  
9012-8774  
9012-8780  
2010-2238STD

9250-1180  
9250-1025  
9012-8774  
9012-8780  
2010-2238CE

\* NOM - Specifications are subject to change without notice