

ASKIR C30

ON CASTORS

A surgical suction unit indicated for professional use.

Thanks to its powerful aspiration and to the double collection system it finds application in general surgery and

The device is suitable for nasal, oral and endotracheal aspiration of bodily fluids (mucus or catarrh) from adults and children.

MAIN FEATURES

Non-stop operation

In order to guarantee the best performance, engine can work without interruption 24/24h.

Maneuvrable

Conveniently placed on a 5-castors stand to make its displacement easier.

High-quality materials

The chassis is made of UL94-V0 durable, flame retardant, corrosion resistant and rust-proof ABS plastic.

Long lasting accessories

All accessories are made in resistant, autoclavable materials to guarantee safety and hygiene over time.

Best Seller configuration

ASKIR C30 +2000 mL reusable jar in Makrolon Polycarbonate (121°C)
Product code: RE 410250

Select your configuration at page 46

The warranty lapses immediately if the unit is used without all the required accessories or with non-original CA-MI accessories.



APPLICATIONS

GENERAL SURGERY	NEUROSURGERY
GINECOLOGY	OBSTETRICS
OPERATING THEATRE	DENTAL PRACTICE

TECHNICAL DATA

MOTOR	Oiless and maintenance-free piston pump
POWER FEEDING	230 V / 50 Hz
ISO 10079-1 CLASSIFICATION	High vacuum / High flow
MAX VACCUM adjustable	-0,80 Bar -80 kPa -600 mmHg value at sea level different altitudes may affect it
FREE AIR FLOW RATE	40 L/min
NOISE LEVEL	61,5 dB(A)
DUTY CICLE	Non-stop operation
WEIGHT	6,20 Kg
DIMENSIONS	32 x 30 x 99 cm
WARRANTY	2 years
SHIPPING CARTON	1 unit
MANUFACTURING COUNTRY	Italy



INCLUDED ACCESSORIES

Basic Setup

2 REUSABLE JARS OR DISPOSABLE SYSTEM
depending on configuration

AUTOCLAVABLE SILICONE TUBES
ø 8/14 mm

ANTIBACTERIAL & HYDROPHOBIC FILTER
single-patient

BI-CONICAL CONNECTOR
ø 10/11/12 mm

POWER CORD WITH SCHUKO PLUG

FS Setup

ALL THE ABOVE-MENTIONED ACCESSORIES
+ FOOT-SWITCH CONTROL
with intermittent or continuous operation

OPTIONAL ACCESSORIES

Reusable Jars in Makrolon®, APEC® Polycarbonate or Polysulfone

AVAILABLE SIZES:
2000 mL

AVAILABLE SETS:
› Complete set (Jar+cover with overflow system+ O-ring)
› Jar (without cover)
› Cover with overflow system
› O-ring

Further details at page: 41




OPERATING WITH:
REUSABLE COMPLETE TUBES' SET
3 silicone tubes (Ø 8/14mm, L: 150cm/ 18cm each)
+ 1 bi-conical connector Ø 10-11-12mm
+ 1 Antibacterial and hydrophobic filter

Disposable Flovac® System

AVAILABLE SIZES:
2000 mL

AVAILABLE SETS:
› Reusable Flovac® Canister
› Disposable Flovac® Liner
› Male Connectors for Flovac® Liner (5's pack)

Further details at page: 43



OPERATING WITH:
FLOVAC COMPLETE TUBES' SET
2 tubes (Ø 8 / 14mm, L: 140cm / L: 26cm)
+ 1 bi-conical connector (Ø 10/11/12mm)

Rolls of silicone tube

AVAILABLE SIZES:
› Ø 8 / 14mm, Length = 1m
› Ø 8 / 14mm, Length = 10m
› Ø 8 / 14mm, Length = 50m

Cannulas

AVAILABLE MODELS

› Yankauer Handle Flat Tip With Hole
› Yankauer Handle Crown Tip With Hole
› Yankauer Tube L= 180cm

Foot-switch control

with plug-in system

