

EL12D

ELECTRIC TILT TABLE



CONFY CARE®

Electric Tilt Table mainly used in the nervous system out of control caused by weakness, pain, spasm of patients, auxiliary training to rehabilitation purposes.



Features

- Electric tilt operation by handle switch
- Tilt angle inclinator included
- Activity tray & patient straps included
- Sleek design with rounded frame legs for safety
- 125mm retractable castors for safety and mobility
- Accessories: adjustable belts (3pcs), adjustable foot plate, active table adjustable
- Reliable TIMOTION motor
- Padding thickness 5cm
- Foam Density: 36kg/m3
- Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant
- Warranty: 5 Years Construction, 3 Years Electrics
- Certification: CE, FDA, ISO 13485:2016
- Weight Capacity: 250KG



1 Motor



Warranty



Certificate



Max Loading

Specification



190cm



64cm



51cm



-2° to +85°



0° to +20°



5cm



φ 125mm



80kgs



1.16



215 x 90 x 60cm