

## MIKO

### HOME CARE COMPRESSOR NEBULIZER



#### Non-stop operation

The device can work without interruption, in order to guarantee the best performance.



#### Fast & efficient

Fast HI-FLO Nebulizer, 7ml capacity, takes only 6 minutes with 2,5 ml of drug (0,35 ml/min).



#### Compact Design

Ideal for easy handling, equipped with carrying bag for storing accessories and medicines.



#### Quiet

Noise level as low as 56 dB(A).



#### Excellent nebulization

MMAD: 2,44  $\mu$ m, respirable fraction >80%, residual volume: 0,40 -0,50 ml.

## INCLUDED ACCESSORIES

#### HI-FLO

jet nebulizer (see page 10)

#### ADULT MASK + CHILD MASK

with elastic strings

#### MOUTHPIECE + NOSEPIECE

non-invasive

#### AIR TUBE

#### AIR FILTER

#### CARRYING BAG

for safe storage and transport

**PRODUCT CODE** RE 300600/03