

# ELECTRA

## Electrotherapy

Electra delivers different types of currents suited for most of the therapy of painful syndromes and for invigorating or hypertrophying treatment through passive gymnastics, and rehabilitation of normoinnervated and/or denervated muscular tissues.

Electra, realized with the more advanced electronic technologies, it is proposed in 2 or 4 channels.

## Features

- Max output direct current 30 mA
- Max output pulse current 70 mA
- Max output voltage 100 Vdc



Diadynamic



Faradic




Galvanic



Impulsive



T.E.N.S.



Bipolar

## +Whats **New**

- Wide range of preset protocols
- Customized programs
- Planning sequence of treatment phases each channel allowed
- Completely independent channels
- Methods of safety control with deactivation of the output current in case of accidental detachment of the electrodes



**FOR  
MORE  
INFO**