

# PRESSA

## Pressotherapy

Pressa is a modern and reliable compressive system of the new generation. The equipment is destined to physiotherapeutic treatments in which is in demand mechanical pressure massage for mechanical lymphatic draining.

The pressa measurement exactness, the polyfuntionality of the planning system and the action of the special physiological accessories, with overlapping sectors, allow to effect the lymph nodes unblocking before make the pressure massage.

## Features

- Max output pressure 250 mmHg
- Sector pressure range 30-90 mmHg
- Number of sectors 8
- Treatment mode Time or Cycles
- Power supply 115 or 230 Vac
- Electrical input power 70 VA

## +Whats **New**

- Stored therapeutic protocols for each accessories
- Planning sequence of treatment phases
- Silent compressor
- Automatic deflating at the treatment end
- Fast accesories connection



**FOR  
MORE  
INFO**