

## SEA-23L/18L



### AIR FILTER

Prevents airborne impurities from entering the chamber to prevent contamination.

### CREATING NEW CHAMBER STRUCTURE

The chamber features advanced structure, manufactured using advanced shaping, fine drawing, dusting, and polishing techniques. It offers exceptional pressure resistance. The outer layer is formed using ultra-fine microporous foaming technology.

### BIG CAPACITY OPEN TANK

- 1) Open tank, easy-to-clean.
- 2) Built-in water quality sensor monitors water quality, reduces machine failure rate, and prolongs machine life time.

### SUITABLE FOR DIFFERENT TYPES INSTRUMENTS

Surgical instruments, dental instruments, glass wares, injection units and other objects of compression resistance and hot steam resistance.

### DOUBLE LOCK SYSTEM -MECHANICAL LOCK WITH AIR LOCK SYSTEM

Double lock system makes the autoclave safer, and prevents potential safety hazards