




Fetal Monitor

BT-350

- 7" wide TFT color LCD display
- Up to 150 patients data saving
- USB data saving
- Multiple language support
- Built-in quick guide [LCD Type only]
- Waterproof probe
- GE / Philips paper compatible
- Interpretation function



LED Type



LCD Type



Fetal Monitor **BT-350**

Specifications

FHR	• Dual pulsed doppler • Intensity : 10mW / cm² or less • Accuracy : $\pm 2\%$ of range	• Ultrasound frequency : 1 MHz • FHR range : 30~240 bpm FHR 50~210 bpm FHR
FM	• Auto-detection of dual fetal movement	
UC	• External type • Reference[zero] control	• Frequency response : DC~0.5HZ • Measurement range : 0~99units
Printer	• Thermal array type • Auto print : off, 10, 20, 30, 40, 50, 60 min	• Print speed : 1, 2, 3cm/min, high speed(Trend mode) • Paper feeding function
Display	• 3 Channels(FHR I , FHR II , UC) • 7" Wide TFT-color LCD : LCD type • Large 7-segment LED type	
Function	• Mark function • Auto print function • Interpretation function	• FHR II offset function • Fetal movement print function • Multi-language & quick guide function : LCD type
Trend	• USB Data Saving	• Data saving for 450 hours(3hr/person) : LCD type only
Power	• Input : AC 100~240V(50/60Hz) • Battery operation(Optional)	• Output : DC 18V, 2.5A • Consumption : 80VA
Dimension	• 80 mm(H) x 330 mm(L) x 280mm(D)	
External link	• Wire : RS-232C	• Wireless : Bluetooth
Accessory standard	• Ultrasound doppler probe 2ea • Event marker jack 1ea • Power adaptor & Cord 1ea • Probe belt 3ea	• UC probe 1ea • Print paper 2ea • Ultrasound gel 1ea • Operation manual 1ea • Probe holder 2ea
Optional	• AST(Acoustic stimulator) • Bluetooth system	• Rechargeable battery • Wall mount bracket • Carrying bag
Warranty	• 1 year (Accessory excluded)	

PC Interface (Central monitoring system)

BCM-350

- Remote control up to 8 monitors
- Fetal heart sound record & play
- Unlimited multi language support
- Wire(RS-232C)and wireless(Bluetooth) connection
- Network system by LAN connection(BCM-350 only)
- Remote real-time picking-up fetal heart sound(Patent)


Installation



Notebook style



Wall mount style



7th Fl., A Bldg., Woolim Lions Valley 5-cha, 144-3, Sangdaewon-dong, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea
 Tel 82-31-750-0340 Fax 82-31-750-0344
 E-mail bistos@bistos.co.kr <http://www.bistos.co.kr>