

SYRINGE PUMP

SIMPLE AND SAFE



QS-800

- Syringe handle clip makes loading easier to Replace syringe easily and avoid siphon.
- Electric control and Auto-induction.
- Accurate control for smooth flow.
- 0.01 ml/h increment for high risk drug transfusion.
- Extra visible and graphical alarms.



GENERAL SPECIFICATIONS

Syringe sizes	5ml,10ml,20ml,30ml,50/60ml
Flow rate accuracy	± 2%
Flow rate	0.01-1500 ml/h
Bolus settings	Manual and programmable function
KVO rate	0.01-5.00 ml/h
Infusion modes	7 infusion mode
Occlusion levels	12 levels Four kinds of unit selectable (mmHg,kPa,bar,psi)
Alarm	VTBI near end, VTBI infused, Pressure high, Battery nearly empty, Battery empty, No power supply, Reminder alarm, KVO finished
Battery life	Approx.12 hours(at 5ml/h flow rate with new fully charged battery)
Weight	Approx.1.8kg