

Model: 3001B Hydraulic Operation Table (Side operating)



Applications:

The operating table is supplied for thoracic surgery, abdominal surgery, ophthalmology, otorhinolaryngology, obstetrics and gynecology, urinary surgery and orthopaedics.

Feature:

- base shell, stand column and bed surface is made of stainless steel.
- Side operating universal table.
- Lifting of this operating table is controlled by foot oil pump, and the rest of the movements are controlled by mechanical transmission side manipulation.
- Shoulder support, body support, arm rest, with leather wrap.

Specification:

1. Table Length (mm): 2100±50
2. Table Width (mm): 500±20
3. Table Height Range (mm): (750~1000)±50
4. Tilt(Left/Right): ≥18°/≥18°
5. Tilt(forward/backward): ≥30°/≥20°
6. Headboard upwards: ≥70°, Downwards: ≥90°
7. Backboard upwards: ≥75°, Downwards: ≥18°
8. Leg board downwards: ≥90°
9. Waist board lifting: 120±10mm

Standard Accessories:

A pair of shoulder support frame, a pair of body support frame, a pair of arm support frame, a pair of leg support frame, a screen frame