

FL-A Series X-Ray Film Viewer

Feature

The average luminance of viewing screen increase from 1200cd/m² to over 2000cd/m²,
The uniformity is over 70%, higher and more uniformity luminance shows the ample imaging hierarchy,
and provide more diagnostic information for doctors.
High intensity aluminum-alloy of the frame and body, keep improving technics, Firm and durable.

Specification:

1. Fixed luminance, ON/OFF button
2. Average luminance: $\geq 2000 \text{ cd/m}^2$
3. Average illuminance: $\geq 6300 \text{ lux}$
4. Thinness: 10cm
5. Lamp: LED lamp bar



Single connection

Screen Dimension: 360*435mm
Shape Dimension: 440*580*100mm
Consume Power: ≤ 50

Double connection

Screen Dimension: 720*435mm
Shape Dimension: 800*580*120mm
Consume Power: ≤ 80

Triple connection

Screen Dimension: 1080*435mm
Shape Dimension: 1160*580*120mm
Consume Power: ≤ 120