

TR2013-620 Ceiling Type LED Operation Lamp



Features:

1. Adopt imported CREE LED cold light source, high quality, stable and reliable.
2. LED light has the large Illumination, better than ordinary light bulbs, save energy over 75%, economical and practical.
3. Ceiling type operating light is 360° range of circular provides for an easy full range of motion;
4. Advanced dimmer system which allows for easy manipulation to meet the user's needs for illumination;
5. The streamline lamp is designed according to the dynamics of its atmosphere to provide a sterile operating environment;

Technical Parameter:

1. Diameter of lamp head: 670mm
2. Number of LED bulbs: 46pcs
3. Illuminance: $\geq 160000\text{Lux}$
4. Power supply: $\sim 220\text{V}$ 50Hz
5. Input power: 400VA
6. Color indication index: $\geq 97\%$
7. Bulb life: 250000h
8. Six-dimensional joint: 360°
9. Lowest height of installation: 2600mm