

TR2015D-350 OPERATION LIGHT



Technical parameter

Power supply :	~ 220V
LED bulb power:	DC:24V
Luminance:	≥70.000Lux
Input power:	90VA
Color temperature :	4000±500K
Spot diameter:	150~300cm
Environmental temperature:	-40℃ ~ +55℃
Relative humidity:	≤90%
Fuse:	RT1-20 T2A L250V RF1-20 F5A L250V
LED bulb power:	(3W×21pcs) 63W
Color indication index:	≥85
Weight:	29Kg