

TRLEDSZ4 LED Operation Light Single head ceiling type

Feature:

- ◆ Imported LED cold light source, with the ultrathin optical lens.
- ◆ Digital means are adopted for stepless regulation of LED brightness.
- ◆ Light field size can be adjusted by the handle.
- ◆ The removable handle sheath can be disinfected.
- ◆ The service life of LED shadowless lamp is long (50000hours).
- ◆ The no welding titanium alloy arm adopt the imported spring, make the arm lift and lower easily and freely.



Main Technical Parameter.

Model	TRLEDSZ4
Type	Single head ceiling type
Diameter of light head	
Illumination at 1m (lux)	$100000 \geq EC \geq 40000$
Light field diameter at 1m (cm)	10-28
Diameter of spot D_{10} (mm)	180 ± 40
Diameter of spot D_{50} (mm)	110 ± 35
L1+L2 Depth of illumination (cm)	-
Color rendering index (Ra)	$100 \geq Ra \geq 85$
Color temperature (K)	4500 ± 500
Number of LED bulbs	32
Average bulb life (hours)	50000 hours
Heat-to-light ratio ($mW/m^2 \text{ lx}$)	$< 6 mW/m^2 \text{ lx}$
Working distance (cm)	70-140
Power supply	100-240V, 50/60Hz
Input power (VA)	60
Ceiling height for installation (cm)	270-320