

YT-1

Infrared Forehead Thermometer

CE & CFDA APPROVED

• Measurement range :	Human body temperature mode: 32.0°C~43.0°C
• Display Resolution :	0.1°C
• Measurement place :	Forehead
• Accuracy :	32.0°C~34.9°C, maximum permissible error $\pm 0.3^{\circ}\text{C}$ 35.0°C~42.0°C, maximum permissible error $\pm 0.2^{\circ}\text{C}$ 42.1°C~42.9°C, maximum permissible error $\pm 0.3^{\circ}\text{C}$
• Display screen :	LCD screen with 4-digit and special icons
• Operating environment :	10°C~40°C
• Max. relative humidity :	$\leq 85\%$
• Storage environment :	-20°C~+55°C
• Max. relative humidity :	$\leq 85\%$
• Battery:	2 $\times$ 1.5V AAA batteries

