



Medical Electric thermometer YT306

- ▶ Measure scope: 32.0°C ~ 42.9°C (89.6°F ~ 109°F)
- ▶ Accuracy : $\pm 0.1^{\circ}\text{C}$ ($\pm 0.2^{\circ}\text{F}$)
- ▶ Measuring time: about 1 minute (Mouth use)
2-5 minutes (armpit use)
- ▶ Alarm Function: Thermometer beeps after temperature taken when
temperature higher than 37.8°C,
It will beep in different sounds.
- ▶ Unit of Display: 0.1°C
- ▶ Memory Function: Remember the last measuring result.
- ▶ Power: 0.1 mW
- ▶ Battery: 1X1.5V LR41, Ordinary button battery,
useful life about 1000times.
- ▶ Packing: 1pc/plastic cover/color box, 200pcs/carton,
- ▶ Certifications: ISO 9001:2008, ISO 13485:2003, RoHS European CE Mark



JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Manufacturing Address: YunYang Industrial Park, DanYang, Jiangsu Province, P.R. China

Tel: 86-511-8690 0869 86-511-8690 0859 Fax: 86-511-8690 0860

E-mail: yuwellsales@yuyue.com.cn ; sales@yuyue.com.cn

www.yuwell.com



/globalyuwell