Folding Stretcher Td010115 MLF999A

MYLIFE







Aluminum Alloy Foldable Stretcher (Folding In Four Parts)
MLF999-A(TD010 115)

● HIGH-STRENGTH ALUMINIUM POLES AND VINYL COATED NYLON/WATERPROOF MATERIALS.

● FOLDS LENGTHWISE FOR COMPACT STORAGE.

● LIGHT-WEIGHTED, SMALL-SIZED, EASY-CARRIED, USE-SAFETY.

Portable Collapsible Pole Emergency Folding Stretcher (TD010115)

*High-strength Aluminum poles and vinyl coated /waterproof materials.
*Folds lengthwise for compact storage
*Light-weighted, small-sized, easy-carried, use-safety.
*Commonly used in hospital, gymnasium, ambulance, and battle fields to carry patients and wounded person.

Model	N. W	G. W	Load Weight	Dimension of packing box	Product Size	
MLF999A(TD010115)	7kg	8kg	159kg	110*52*22cm(2pcs/box)	190*48*20cm	